



06/05/2025

K. Varunkumar
BE23B016

NOTE: The above schematic isn't to scale; proper dimension are provided with the fabrication process; a zoomed version is used in the process diagram below

i) Take the Si wafer, clean it using Piranha soln, followed by clean rinse in DI water



ii) SiO_2 formation, using thermal oxidation of approx 2 μm



iii) a) primer treatment (HMDS) ; b) Spin coating of -ve PR ;
c) Soft bake ; d) Align mask shown in fig (c) with the setup ;

e) UV ~~treatment~~ exposure ; f) PR development followed by hard bake in a hot plate ; g) Standard cleaning ; h) Deposition of ITO via E-beam evaporation ; i) Removal of PR (Acetone) → Lift-off

Now the ITO electrodes have been patterned.




iv) Primer treatment → SU-8 spin coating (HMDS)

Mask Alignment [Figure (e)] ← Soft bake

UV exposure → PR development

Standard clean. ← Hard bake




Now the cell traps are made using SU-8; → figure (f)

v) The substrate is now flipped over to ensure the through can be drilled.

Primer treatment → SU-8 spin coating (HMDS)


Mask Alignment Figure (g) ← Soft bake



(TMR, diameter = 10 μm)

The section is through the midpoint/center of the cell trap, thus ITO and SU-8 aren't visible

Contd..



NOTE: PR is removed using tweezers to ensure the SU-8 trap stays intact / Alternatively, only the bottom PR can be exposed to Acetone for its selective removal.

vi) The top is sealed with a glass that reduces the headspace for the fluid flow through the narrow channels as shown in fig (h).

The bottom is sealed with glass, which contains Titanium micro-rings for photoporation.

Working mechanism :

- This device integrates, photo-poration with electroporation to ensure high transfection efficiency, as well as high cell viability.

- Cells are trapped in nearly spherical regions of 35 μm dia, using hydrodynamics, the gap channels ensure low resistance when uncovered & high resistance when covered,

This along with the appropriate dimensioning ensure high efficiency single cell trapping.

- Once the cells ~~are~~ are positioned, from under the substrate IR light of approx 1050 nm is pulsed at regular intervals. (The wavelength also ensure minimal cell damage $E \propto 1/\lambda$)

The laser excites the TiO_2 ring which concentrates the energy and forms photothermal bubbles. These photothermal bubbles then lyse the cell wall, which allows drugs to enter via diffusion.

- Finally we have the ITO nanogap electrodes. The nanogap, ensure a concentrated electric field near the tip of the electrode, ~~cause~~ causing the cell-wall lysis in it's locality. Drugs once again enter the cell through this via diffusion.

Expected outcomes

- Certain methods of drug delivery prove to be inefficient against certain cell types, however, by combining 2 techniques, we can ensure regardless of cell type, there is a high transfection efficiency

- Both methods ~~are~~ combined here are very localized and as such don't lead to multiple pores generating simultaneously which can cause cell death.

- Through the combination of the 2 techniques, if a particular method is found to be more effective, it can carry the larger burden. This compact design also allows for 1 million cells in less than 100 cm^2 area.